


## Product Datasheet

### Human SERPINA3 ELISA Kit (orb562027)

<b>Description</b>	Human SERPINA3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100nmol/ml
<b>Target</b>	SERPINA3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P01011</b>
<b>Sensitivity</b>	0.938nmol/ml
<b>Expiration Date</b>	Please enquire.



Human  
SERPINA3(Alpha-1-  
antichymotrypsin)...