


## Product Datasheet

### Human PLG ELISA Kit (orb562038)

<b>Description</b>	Human PLG ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.125-5000pg/ml
<b>Target</b>	Plg
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P00747</b>
<b>Sensitivity</b>	46.875pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PLG(Plasminogen)  
ELISA Kit