


## Product Datasheet

### Human Caspase 1 ELISA Kit (orb562045)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Human Caspase 1 ELISA Kit  |
| <b>Reactivity</b>      | Human  |
| <b>Range</b>           | 78.125-5000pg/ml   |
| <b>Target</b>          | CASP1  |
| <b>Storage</b>         | 2-8 °C for 6 months  |
| <b>Note</b>            | For research use only  |
| <b>Sample Types</b>    | Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples |
| <b>Uniprot ID</b>      | <b>P29466</b>  |
| <b>Sensitivity</b>     | 46.875pg/ml  |
| <b>Expiration Date</b> | Please enquire.  |



Human  
CASP1(Caspase-  
1) ELISA Kit