


## Product Datasheet

### Human PDGFD ELISA Kit (orb562055)

<b>Description</b>	Human PDGFD ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	12.5-800pg/ml
<b>Target</b>	PDGFD
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9GZP0</b>
<b>Sensitivity</b>	7.5pg/ml
<b>Expiration Date</b>	Please enquire.



Human PDGFD(Platelet-derived growth fact...