


Product Datasheet

Human IL36G ELISA Kit (orb562065)

Description	Human IL36G ELISA Kit
Reactivity	Human
Range	31.2-2000pg/ml
Target	IL36G
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9NZH8
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
IL36G(Interleukin-
36 gamma) ELISA
...