


## Product Datasheet

### Human CCL28 ELISA Kit (orb562097)

<b>Description</b>	Human CCL28 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	CCL28
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9NRJ3</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
CCL28(C-C  
motif  
chemokine  
28) ELIS...