


## Product Datasheet

### Human VIP ELISA Kit (orb562108)

<b>Description</b>	Human VIP ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.313-20pg/ml
<b>Target</b>	VIP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P01282</b>
<b>Sensitivity</b>	0.188pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
VIP(VIP  
peptides)  
ELISA Kit