


Product Datasheet

Human SERPINA6 ELISA Kit (orb562155)

| | |
|------------------------|--------------------------|
| Description | Human SERPINA6 ELISA Kit |
| Reactivity | Human |
| Range | 15.625-1000ng/ml |
| Target | SERPINA6 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P08185 |
| Sensitivity | 9.375ng/ml |
| Expiration Date | Please enquire. |



Human
SERPINA6(Corticosteroid-
binding gl...