


## Product Datasheet

### Human Adrenomedullin ELISA Kit (orb562176)

<b>Description</b>	Human ADM(Adrenomedullin) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	6.25-400pg/ml
<b>Target</b>	ADM
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P35318</b>
<b>Sensitivity</b>	3.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
ADM(ADM)  
ELISA Kit