


## Product Datasheet

### Human Enkephalin ELISA Kit (orb562209)

<b>Description</b>	Human Enkephalin ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78-5000pg/ml
<b>Target</b>	PENK
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P01210</b>
<b>Sensitivity</b>	46.9pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
PENK(Proenkephalin-  
A) ELISA Kit