


Product Datasheet

Human DEFB4A ELISA Kit (orb562226)

Description	Human DEFB4A ELISA Kit
Reactivity	Human
Range	62.5-4000pg/ml
Target	DEFB4A
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	015263
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Human
DEFB4A(Beta-
defensin 4A)
ELISA Kit