


## Product Datasheet

### Human NCOR2 ELISA Kit (orb562244)

<b>Description</b>	Human NCOR2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.125-8ng/ml
<b>Target</b>	NCOR2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9Y618</b>
<b>Sensitivity</b>	0.075ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
NCOR2(Nuclear  
receptor  
corepressor...