

Product Datasheet

Human TAC1(Protachykinin-1) ELISA Kit (orb562267)

Description Human TAC1(Protachykinin-1) ELISA Kit

Reactivity Human

Range 4.69-300pg/ml

Target TAC1

Storage 2-8 °C for 6 months


Note For research use only

Application notes TAC1 ELISA Kit allows for the in vitro quantitative determination of TAC1 concentrations in serum, plasma, tissue homogenates and other biological fluids.

Uniprot ID **P20366**

Sensitivity 2.813pg/ml

Expiration Date Please enquire.



Human TAC1(Protachykinin-1) ELISA Kit