


## Product Datasheet

### Human SHBG ELISA Kit (orb562268)

<b>Description</b>	Human SHBG ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pmol/L
<b>Target</b>	SHBG
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P04278</b>
<b>Sensitivity</b>	9.375pmol/L
<b>Expiration Date</b>	Please enquire.



Human SHBG(Sex hormone-binding globulin)...