


## Product Datasheet

### Human PDE5A ELISA Kit (orb562288)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Human PDE5A ELISA Kit  |
| <b>Reactivity</b>      | Human                  |
| <b>Range</b>           | 0.156-10ng/ml          |
| <b>Target</b>          | PDE5A                  |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">O76074</a> |
| <b>Sensitivity</b>     | 0.094ng/ml             |
| <b>Expiration Date</b> | Please enquire.        |



Human  
PDE5A(cGMP-  
specific 3',  
5'-cyclic ...