


## Product Datasheet

### Human SLC25A5 ELISA Kit (orb562307)

<b>Description</b>	Human SLC25A5 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78-5000pg/ml
<b>Target</b>	SLC25A5
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P05141</b>
<b>Sensitivity</b>	46.9pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
SLC25A5(ADP/ATP  
translocase 2)  
ELI...