


Product Datasheet

Human UBQLN2 ELISA Kit (orb562313)

Description	Human UBQLN2 ELISA Kit
Reactivity	Human
Range	78-5000pg/ml
Target	UBQLN2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9UHD9
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Human
UBQLN2(Ubiquilin-
2) ELISA Kit