


## Product Datasheet

### Human Smad1 ELISA Kit (orb562361)

<b>Description</b>	Human Smad1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.781-50ng/ml
<b>Target</b>	Smad1/MADH1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q15797</a>
<b>Sensitivity</b>	0.469ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
Smad1/MADH1(Mothers  
Against Decape...