


## Product Datasheet

### Human KLK8 ELISA Kit (orb562378)

<b>Description</b>	Human KLK8 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	KLK8
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>O60259</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
KLK8(Kallikrein-  
8) ELISA Kit