


## Product Datasheet

### Human CA2 ELISA Kit (orb562414)

<b>Description</b>	Human CA2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.125-200ng/ml
<b>Target</b>	CA2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P00918</b>
<b>Sensitivity</b>	1.875ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
CA2(Carbonic  
anhydrase 2)  
ELISA Ki...