


## Product Datasheet

### Human IMA ELISA Kit (orb562428)

<b>Description</b>	Human IMA ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	7.8-500ng/ml
<b>Target</b>	IMA
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	4.688ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
IMA(Ischemia  
Modified  
Albumin)  
ELI...