


## Product Datasheet

### Human PDIA4 ELISA Kit (orb562576)

<b>Description</b>	Human PDIA4 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/ml
<b>Target</b>	PDIA4
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P13667</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Human PDIA4(Protein disulfide-isomerase ...