


## Product Datasheet

### Human AFP-L3 ELISA Kit (orb562664)

<b>Description</b>	Human AFP-L3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	AFP-L3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
AFP-  
L3(Alpha-  
fetoprotein  
Lens culi...