


Product Datasheet

Human SYNPO ELISA Kit (orb562676)

Description	Human SYNPO ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	SYNPO
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Cell or tissue lysate,Other liquid samples
Uniprot ID	Q8N3V7
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
SYNPO(Synaptopodin)
ELISA Kit