


## Product Datasheet

### Human CNDP1 ELISA Kit (orb562774)

<b>Description</b>	Human CNDP1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/ml
<b>Target</b>	CNDP1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<b>Q96KN2</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Human CNDP1(Beta-Ala-His dipeptidase) EL...