


## Product Datasheet

### Human NTs ELISA Kit (orb562796)

<b>Description</b>	Human NTs ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100pg/ml
<b>Target</b>	NTs
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Uniprot ID</b>	<b>P30990</b>
<b>Sensitivity</b>	0.938pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
NTs(Neurotensin)  
ELISA Kit