


Product Datasheet

Human MTRNR2L6 ELISA Kit (orb562859)

Description	Human MTRNR2L6 ELISA Kit
Reactivity	Human
Range	0.312-20ng/ml
Target	MTRNR2L6
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P0CJ73
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
MTRNR2L6(Humanin-
like protein 6) E...