


Product Datasheet

Human PIGR ELISA Kit (orb562888)

| | |
|------------------------|-----------------------|
| Description | Human PIGR ELISA Kit |
| Reactivity | Human |
| Range | 31.25-2000pg/ml |
| Target | PIGR |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P01833 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Human
PIGR(Polymeric
Immunoglobulin
Rece...