


Product Datasheet

Human YWHAQ ELISA Kit (orb563044)

Description	Human YWHAQ ELISA Kit
Reactivity	Human
Range	0.78-50ng/ml
Target	YWHAQ
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	P27348
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Human
YWHAQ(14-
3-3 protein
theta) ELISA
...