


## Product Datasheet

### Human BTG3 ELISA Kit (orb563076)

<b>Description</b>	Human BTG3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.2-2000pg/ml
<b>Target</b>	BTG3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q14201</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
BTG3(Protein  
BTG3) ELISA  
Kit