


Product Datasheet

Human DDR1 ELISA Kit (orb563206)

| | |
|------------------------|------------------------|
| Description | Human DDR1 ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/ml |
| Target | DDR1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q08345 |
| Sensitivity | 0.094ng/ml |
| Expiration Date | Please enquire. |



Human
DDR1(Epithelial
discoidin
domain-c...