


## Product Datasheet

### Human MPEG1 ELISA Kit (orb563213)

<b>Description</b>	Human MPEG1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/ml
<b>Target</b>	MPEG1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich ELISA, Double Antibody
<b>Uniprot ID</b>	<b>Q2M385</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
MPEG1(Macrophage-  
expressed gene 1 ...