


## Product Datasheet

### Human HSPA5 ELISA Kit (orb563238)

<b>Description</b>	Human HSPA5 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625-40ng/ml
<b>Target</b>	HSPA5
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P11021</b>
<b>Sensitivity</b>	0.375ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
HSPA5(78  
kDa  
glucose-  
regulated  
pro...