


## Product Datasheet

### Human ARHGAP21 ELISA Kit (orb563257)

<b>Description</b>	Human ARHGAP21 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/ml
<b>Target</b>	ARHGAP21
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q5T5U3</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
ARHGAP21(Rho  
GTPase-  
activating  
pro...