


## Product Datasheet

### Human MEP1B ELISA Kit (orb563282)

<b>Description</b>	Human MEP1B ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/ml
<b>Target</b>	MEP1B
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q16820</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
MEP1B(Meprin  
A subunit  
beta) ELISA...