


## Product Datasheet

### Human SFTPB ELISA Kit (orb563284)

<b>Description</b>	Human SFTPB ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.781-50ng/ml
<b>Target</b>	SFTPB
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P07988</b>
<b>Sensitivity</b>	0.469ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
SFTPB(Pulmonary  
surfactant-  
associa...