


## Product Datasheet

### Human PRAP1 ELISA Kit (orb563287)

<b>Description</b>	Human PRAP1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-200pg/ml
<b>Target</b>	PRAP1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q96NZ9</b>
<b>Sensitivity</b>	1.875pg/ml
<b>Expiration Date</b>	Please enquire.



Human PRAP1(Proline-rich acidic protein ...