


## Product Datasheet

### Human CTXIII ELISA Kit (orb563304)

<b>Description</b>	Human CTXIII(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	CTXIII
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human CTXIII(C-terminal telopeptide of C...