

Product Datasheet

FANCG (phospho-Ser383) antibody (orb5634)

Description	FANCG (phospho-Ser383) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FANCG around the phosphorylation site of Ser383 (RF(p-S)PP)
Target	phospho-FANCG (Ser383)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69 kDa
Uniprot ID	O15287
Expiration Date	12 months from date of receipt.