


Product Datasheet

Human ANGPTL2 ELISA Kit (orb563409)

Description	Human ANGPTL2 ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	ANGPTL2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q9UKU9
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
ANGPTL2(Angiopoietin-
related prote...