


## Product Datasheet

### Human CCR6 ELISA Kit (orb563451)

<b>Description</b>	Human CCR6 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/ml
<b>Target</b>	CCR6
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P51684</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
CCR6(C-C  
chemokine  
receptor  
type 6...