


Product Datasheet

Human 4-HNE ELISA Kit (orb563477)

Description	Human 4-HNE ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	4-HNE
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,urine,Other liquid samples
Sensitivity	18.75 pg/ml
Expiration Date	Please enquire.



Human 4-HNE(4-Hydroxynonenal) ELISA Kit