


Product Datasheet

Human MCL1 ELISA Kit (orb563493)

| | |
|------------------------|-----------------------|
| Description | Human MCL1 ELISA Kit |
| Reactivity | Human |
| Range | 31.2-2000pg/ml |
| Target | MCL1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q07820 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Human
MCL1(Induced
myeloid
leukemia
cell...