


## Product Datasheet

### Human SIRT6 ELISA Kit (orb563510)

<b>Description</b>	Human SIRT6 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78-5000pg/ml
<b>Target</b>	SIRT6
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q8N6T7</a>
<b>Sensitivity</b>	46.9pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
SIRT6(NAD-  
dependent  
deacetylase  
si...