


Product Datasheet

Human MEP1A ELISA Kit (orb563559)

Description	Human MEP1A ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	MEP1A
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
Uniprot ID	Q16819
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
MEP1A(Meprin
A subunit
alpha) ELIS...