


Product Datasheet

Human SMOX ELISA Kit (orb563578)

| | |
|------------------------|-----------------------|
| Description | Human SMOX ELISA Kit |
| Reactivity | Human |
| Range | 0.78-50ng/ml |
| Target | SMOX |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q9NWM0 |
| Sensitivity | 0.469ng/ml |
| Expiration Date | Please enquire. |



Human
SMOX(Spermine
oxidase) ELISA
Kit