


Product Datasheet

Human FOLR2 ELISA Kit (orb563613)

| | |
|------------------------|-----------------------|
| Description | Human FOLR2 ELISA Kit |
| Reactivity | Human |
| Range | 15.6-1000pg/ml |
| Target | FOLR2 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P14207 |
| Sensitivity | 9.375pg/ml |
| Expiration Date | Please enquire. |



Human FOLR2(Folate receptor beta) ELISA

...