


## Product Datasheet

### Human CD5L ELISA Kit (orb563617)

<b>Description</b>	Human CD5L ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	CD5L
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>O43866</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
CD5L(CD5  
antigen-  
like)  
ELISA Kit