


## Product Datasheet

### Human CFHR1 ELISA Kit (orb563618)

<b>Description</b>	Human CFHR1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.78-50ng/ml
<b>Target</b>	CFHR1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q03591</a>
<b>Sensitivity</b>	0.469ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
CFHR1(Complement  
factor H-related ...