


## Product Datasheet

### Human GFRA1 ELISA Kit (orb563624)

<b>Description</b>	Human GFRA1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	GFRA1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P56159</b>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human  
GFRA1(GDNF  
family  
receptor  
alpha-1...